



## e-chamber Spacer range

### PHARMAC Listed Respiratory Devices

- Anti-Static ● Reusable, single patient use ● Latex & BPA free
- Made with environmentally friendly durable materials
- Unique patented butterfly one-way valve
- Universal 15mm - 22mm mouthpiece
- Universal metered dose inhaler (MDI)
- Detachable parts - for easy cleaning

BLISTER  
PACKED



ANTI  
STATIC



NO  
PRIMING  
REQUIRED

**Turbo 220ml**  
Neonate, Infant, Child and Adult

**La Grande 510ml**  
Adult

## e-chamber Mask

**S** Small (Infant, Child)

- Contoured
- Anti-collapsible
- Latex & BPA free
- Made from durable and soft thermo-grade silicone
- Compact storage with e-chamber™ spacer range
- Universal 22mm round fit



## BIRD HEALTHCARE

### Introducing the e-chamber™

The **e-chamber™** range of spacers (or holding chambers) with **Butterfly Valve Technology™** are designed to help improve delivery of pressurized metered dose inhaler (puffer) medications into the lungs.

**Proudly distributed by Apex Medical NZ Ltd**



**apexmedical**

P: 0800 469 633  
E: [info@apexmedical.co.nz](mailto:info@apexmedical.co.nz)  
[www.apexmedical.co.nz](http://www.apexmedical.co.nz)



As listed in Section B and Part III of Section H of the Pharmaceutical Schedule, with effect 1 November 2015

Item code	Pharmacode	Description	Order Quantity	Price (GST exclusive)
BHC-TURBO-NZ	2490072	e-chamber Turbo Spacer 220ML	Box/50	\$2.95 ea
BHC-LA-GRANDE-NZ	2490080	e-chamber La Grande Spacer 510ml	Box/25	\$5.12 ea
BHC-MASM-SMALL-NZ	2490099	e-chamber Mask Small	Box/50	\$2.20 ea

To order, please email [spacerorders@apexmedical.co.nz](mailto:spacerorders@apexmedical.co.nz) or phone 0800 273963 (0800 APEXMED)  
Individual consumer sales are available online at Asthma New Zealand [www.asthma-nz-shop.org.nz/Spacers](http://www.asthma-nz-shop.org.nz/Spacers)



apexmedical

[www.apexmedical.co.nz](http://www.apexmedical.co.nz)