

MEDIKRO nano

PC-based Mobile Spirometer

- ▶ Ultimate Mobility
- ▶ Incomparable Durability
- ▶ Patient Safety
- ▶ Environmentally Friendly

ASK
ABOUT OUR
INTRODUCTORY
PRICES



Product Code: M9487

Medikro® Nano is the world's smallest Windows-compatible spirometer, weighing only 11 g. With Nano, the advantages of PC-based spirometry are now available in a mobile setting.

Nano utilizes the easy-to-use Medikro® Spirometry Software (MSS). The system features fulfil the needs of doctors who want to perform on-site testing with a mobile system, while maintaining precision and accuracy of results. With disposable Medikro® SpiroSafe flow transducers there is no need for bacterial filters or sterilization. This increases patient throughput and saves operational costs.



Nano Highlights



Mobility

Medikro Nano is the world's smallest Windows-compatible spirometer, weighing only 11 g. The advantages of PC-based spirometry are now available in a mobile setting.

Small and durable

Minimized electronics reduce hardware failures, and the need for excessive maintenance procedures.



User-friendly

Medikro Spirometry Software is user-friendly and provides intuitive navigation and use.



Patient safety

Disposable SpiroSafe flow transducers provide security in infection control. No exhaled air gets into the pressure tube or the device. This prevents patient-to-patient cross-contamination and protects the spirometer.



Multiple calculation criteria and interpretation options

MSS provides multiple options for calculation of Best PEF, FEV% and FEF. In addition, interpretation options are available for ventilation function, spirometry findings and the bronchodilatation test.

PDF export

Print final reports in PDF format or as images, or let the software do it for you automatically.



Consumables & Accessories



MEDIKRO® SpiroSafe

Product code: M9256



MEDIKRO® Calibration Syringe

Product code: M9474 (3L) / M9477 (1L)



MEDIKRO® Ambi

Product code: M9489



🏠 Kellolahdentie 27, FI-70460 Kuopio, Finland
✉ P.O. Box 54, FI-70101, Kuopio, Finland
☎ +358 44 777 7002
📠 +358 17 283 3300
@ sales@medikro.com
🌐 www.medikro.com

Your local Medikro representative:



P: 0800 469 633
E: info@apexmedical.co.nz
www.apexmedical.co.nz